

Incubator IN and IF with SingleDISPLAY Incubator INplus and IFplus with TwinDISPLAY Natural convection or forced air circulation AtmoCONTROL software

Model sizes: 30 / 55 / 75 / 110 / 160 / 260 / 450 / 750 +20 °C to +80 °C

INCUBATOR I Memmert incubators I are at home in the world of research, pharmaceutics and food technology. The heating and control system are especially optimised for low temperatures of up to +80 °C. To prevent temperature overshoots, temperature is increased within a very narrow control range and kept exactly at the setpoint value. As required, the models with natural convection or with forced air circulation are available.





As little air circulation as possible in the incubator


Forced air convection may destroy the protective layer from moist air that is generated during incubation over the samples. This would lead to dehydration of the culture. In a Memmert incubator, the perfect combination of all-round surface heating and temperature control system ensures that incubation generally takes place without forced air circulation. Provided the chamber is fully loaded and forced air circulation is required, it can be precisely adjusted in 10 % steps from 0 to 100 %.

Sterilisation

The chamber of the incubators INplus/IFplus, including all installations and sensors, can be sterilised at +160 °C in a 4-hour programme to guarantee optimum hygiene.

Fresh air is preheated

Temperature deviations caused by fresh air can influence sample characteristics or prolong drying. In Memmert incubators, the fresh air is therefore fed through a pre-heating chamber and seamlessly introduced into the working chamber.



Air supply from outside

INCUBATORS I

according to DIN 12880:2007-05, EN 61010-1 (IEC 61010-1), EN 61010-2-010

Standard units are safety-approved and bear the test marks: CE







Interior:

Stainless steel, material 1.4301 (ASTM 304) with all-round deep-drawn ribs to integrate the largearea heating with ceramic-metal sheath

Housing:

Textured stainless steel, rear zinc-plated steel, intuitively operated SingleDISPLAY or TwinDISPLAY (TFT colour display) with touchscreen; inner glass door, outside fully insulated stainless steel door (from size 450 two leaves)

Admixture of pre-heated fresh air by electronically adjustable airflap Fresh air:


Connection: Mains cable with plug (German type)

Installation: 4 feet; sizes 450/750 mounted on lockable castors

Interfaces:



USB: only TwinDISPLAY



Model sizes/Descrip	otion		30	55	75	110	160	260	450	750		
Stainless steel	Volume	approx. l	32	53	74	108	161	256	449	749		
interior	Width	(A) mm		400 560 640						1040		
	Height	(B) mm	320	400	560	480	720	800	720	1200		
	Depth (less 39 mm for fan)	(C) mm	250	33	30	40	00	500	60	00		
	Max. number of grids/shelves	number	3	4	6	5	8	9	8	14		
	Max. loading per grid/shelf	kg	20 30							0		
	Max. loading of chamber	kg	60	80	120	175	210		300			
	Max. loading per slide-in drip tray	kg		1,5			3	4		3		
Max. loading per bottom drip tray		kg	1,5				3	4	8			
Textured stainless	Width	(D) mm		585		74	45	824	12	24		
steel exterior	Height (size 450, 750 with castors)	(E) mm	704	784	944	864	1104	1183	1247	1720		
	Depth (without door handle, depth of handle +56 mm)	(F) mm	434	434 514			84	684	684 784			
Standard	Stainless steel grids, electropolished	number	1 2									
equipment	Standard works calibration certificate (measuring point chamber center)	°C	+37									
Temperature	Working temperature range	°C	at least 5 (IN/INplus) 10 (IF/IFplus) above ambient temperature up to +80									
	Setting temperature range	°C	+20 to +80									
	Setting accuracy	°C				0).1					
Further data	Electrical load at 230 V, 50/60 Hz	approx. W	1600	1000	1250	1400	1600	1700	1800	2000		
	Electrical load at 115 V, 50/60 Hz	approx. W	800 900 1500						1500	1800		
Packing data	Net weight	approx. kg	48	57	66	76	96	110	161	217		
J	Gross weight (packed in carton)	approx. kg	64	76	85	101	122	161	227	288		
	Width	approx. mm	660 730 830 930				13	1330				
	Height	approx. mm	890	950	1130	1050	1300	1380	1440	1910		
	Depth	approx. mm	650	6	70	80	00	930	10	50		
Order No. Incubators			IN30	IN55	IN75	IN110	IN160	IN260	IN450	IN750		
I = Incubator			IN30plus	IN55plus	IN75plus	IN110plus	IN160plus	IN260plus	IN450plus	IN750plu		
N = Natural con	vection		IF30	IF55	IF75	IF110	IF160	IF260	IF450	IF750		
F = Forced conv plus = Model with			IF30plus	IF55plus	IF75plus	IF110plus	IF160plus	IF260plus	IF450plus	IF750plu		

Options	30	55	75	110	160	260	450	750	
Voltage 115 V, 50/60 Hz				X2					
Extended overtemperature protection by additionally integrated Pt100 sensor for independent temperature monitoring for models with SingleDISPLAY				A6					
Chamber modification for the application of reinforced perforated stainless steel shelves or stainless steel grids (bearing rails mounted in the working chamber) - includes replacement of standard grids by reinforced grids			-				K	1	
Interior lighting for observing the load				RO					
Interior socket (can only be ordered with limited temperature range - max. +70 °C), ampacity 230 V, 2.2 A, can be switched off with the On/Off switch, cannot be switched individually, moisture tight IP68 (only with SingleDISPLAY), (option A8 necessary)				R3					
Entry port, 23 mm clear diameter, for introducing connections at the side, can be closed by flap and silicone stopper, standard positions				F0 F1 F2 F3					
right centre/top									
Entry port, 23 mm clear diameter, left can be closed by flap, in special				F4					
positions (please state location) right	F5 F6								
Entry port, 14 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location)		D6							
Entry port, 38 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location)	F7								
Entry port, 57 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location; not possible for IF/IFplus size 30-75)		F8							
Entry port, 100 mm clear diameter, can be closed by flap, in special positions in the back wall (please state location; not possible for IF/IFplus size 75)		- F9							
4 - 20 mA current loop interface (0 to +90 °C = 4 - 20 mA) Temperature controller, actual value Temperature of a Pt100 sensor positioned flexibly in chamber for external temperature monitoring (max. 1 SingleDISPLAY, max. 3 TwinDISPLAY) - price per sensor				V3 V6					
Fan speed monitoring with switching-off the heating and with alarm in case of failure - only for IFplus	V4								
Works calibration certificate for 3 temperatures: +37, +52, +70 °C	D00126								
Works calibration certificate for one (freely selectable) temperature value according to customer specification	D00109								
Door with lock and key (safety lock)	B6								
Door hinged on the left		B8 -							
Potential-free contact for combination error message (e.g. supply failure, sensor fault, fuse)	H6								
Door-open-recognition (only for units with TwinDISPLAY)		V5							
MobileALERT, notification by SMS in case of any error or alarm of the device (requires option H6)		G							
Temperature restriction; Temperatures: +50, +55, +60, +65, +70 or +75 °C (Please indicate upon ordering)	A8								
Fresh-air filter (filtration efficiency 80 %) mounted at the appliance bottom (for IF/IFplus, for sizes 30 – 260 requires castor frame R9 or subframe)	R8								
Castor frame (2-part), height 140 mm			R!	9				-	

Accessories	30	55	75	110	160	260	450	750
Stainless steel grid, electropolished		E20164		E20165		E28891	E20182	
Reinforced stainless steel grid, electropolished, max. loading 60 kg; from size 450 with guide bars and fixing screws (only in connection with option K1). Please consider max. loading of chamber	-			E29767		E29766	B32190	
Perforated stainless steel shelf	B29727 B03916 B00325		325	B29725	B00328			
Reinforced stainless steel shelf, max. loading 60 kg; with guide bars and fixing screws (requires option K1). Please consider max. loading of chamber		B32191						
Stainless steel slide-in drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1)	E02070 E02072 E02073		073	E29726	E02075			
Stainless steel slide-in drip tray, 15 mm rim, with guide bars and fixing screws (may affect the temperature distribution, only in connection with option K1)	-						B32	763
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, not in connection with option K1)	B04356 B04358 B04359 B2		B29722	B04	362			
Stainless steel bottom drip tray, 15 mm rim (may affect the temperature distribution, only in connection with option K1)		-					B34055	
Wall bracket for wall mounting		B29756	B29757	B29758	B29759		-	
Guarantee extension by 1 year		GA1Q5						
USB-Ethernet adapter		E06192						
Ethernet connection cable 5 m for computer interface				E06	5189			

Accessories	30	55	75	110	160	260	450	750	
USB User-ID stick (with User-ID licence): Oven-linked authorisation licence (User-ID-programme) on Memory-stick, prevents undesired manipulation by unauthorised third parties. When reordering please specify serial number (only for units with TwinDISPLAY)		B33170							
USB stick with documentation software AtmoCONTROL and operation manual for products with SingleDISPLAY (the standard equipment of appliances with TwinDISPLAY includes one USB stick with AtmoCONTROL). When reordering please specify serial number	B33172								
Set of height adjustable feet (4 pcs)	B29768						-		
Stacking set (4 pcs) for stacking of appliances of same size	B29744						-		
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), straight, for exhaust air ducting (if necessary for connection by hose)	B29718								
Plug-in tube extension (outer diam. 60.3 mm, inner 57 mm), angled, for exhaust air ducting (if necessary for connection by hose)	B29719								
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), with air slots	B29728	B29730	B29732	B29734	B29736	B29738	B29740	B29742	
Flush-fit unit (stainless steel frame covering gap between oven and wall opening), without air slots	B29729 B29731 B29733		B29735	B29737	B29739	B29741	B29743		
Subframe, adjustable in height (size 30 to 75: height 600 mm, size 110 to 450: height 500 mm)	B29745	B29747		B29749		B29751	B29753	-	
Subframe, on castors (size 30 to 75: height 660 mm, size 110 to 160: height 560 mm)	B29746	B29	748	B29	750		-		
Subframe, adjustable in height, height 130 mm, for example for units with fresh air filter		B33	659	B33	661	B33664		-	
Software conforming to FDA AtmoCONTROL. Meets the requirements for the use of electronically stored data sets and electronic signatures as laid down in Regulation 21 CFR Part 11 of the US Food and Drug Administration (FDA). Base licence for the control of one unit (only for units with TwinDISPLAY). Respective IQ/OQ documents available in German and English language (without surcharge)	FDAQ1								
Integration of one additional unit (up to max. 31 units) into an already existent FDA- software licence (only for units with TwinDISPLAY)		FDAQ2							
DAkkS calibration for one free-selectable temperature value according to method C (DKD-R 5-7)		E39696							
DAkkS calibration for further temperature values according to method C (DKD-R 5-7)		E39697							
IQ document with device-specific works test data, OQ/PQ check list as support for validation by customer		D00124							
IQ/OQ document with device-specific works test data for one free-selectable temperature value, incl. temperature distribution survey at Memmert for 9 measuring points (size 30), 27 measuring points (sizes 55 - 1060) to DIN 12880:2007-05. PQ check list as support for validation by customer. 305 € for further temperature values	D00125 D00127								
On-site IQ/OQ for a freely selectable temperature value, including temperature distribution survey for 9 measuring points (size 30), 27 measuring points (sizes 55 - 750) to DIN 12880: 2007-05 (excluding travel costs, not subject to discount, GER, AT, FR only)		DLQ100							
Extension of DLQ100 by an additional freely selectable temperature value (not subject to discount)		DLQ100A							
Individual on-site Performance Qualification (PQ)		DLQ200							
Maintenance UIS - carrying out and documentation according to Memmert maintenance plan (excluding travel costs, not subject to discount, GER, AT, FR only)	S00311								
Maintenance contract UIS - carrying out and documentation according to Memmert maintenance plan, minimum duration 3 years (excluding travel costs, not subject to discount, GER, AT, FR only)				500	311J				
Calibration of one freely selectable temperature value (excluding travel costs, not subject to discount, GER, AT, FR only)	S00205								
Calibration of an additional temperature value (not subject to discount)				S00)215				

PERSONAL NOTES

MODEL VARIANTS

SingleDISPLAY ControlCOCKPIT with one TFT display	TwinDISPLAY ControlCOCKPIT with two TFT displays						
AVAILABLE APPLIANCES	AVAILABLE APPLIANCES						
UN / UF / IN / IF / IPPeco / IPP / UNm / UFm / INm / IFm / SN / SF / IFbw	HPPeco / ICHeco / ICH / HCP / UNplus / UFplus / UF TS / UNpa / VO / INplus / IFplus / ICO / IPPecoplus / IPPplus / ICPeco / ICP / UNmplus / UFmplus / INmplus / IFmplus / SNplus / SFplus / ICOmed						
One high-resolution TFT colour display with touch-sensitive buttons for selection of functions	Two high-resolution TFT colour displays with touch-sensitive buttons for selection of functions						
Available parameters on the ControlCOCKPIT: Temperature (Celsius or Fahrenheit), fan speed, exhaust air flap position, programme time	Available parameters on the Control COCKPIT: All parameters of the SingleDISPLAY and device-specific parameters like relative humidity, illumination and CO ₃						
One temperature sensor Pt100 DIN class A in a 4-wire circuit	Two Pt100 sensors DIN class A in a 4-wire circuit for mutual monitoring, taking over functions in case of an error						
	HeatBALANCE function for application specific adjustment of heat output distribution (balance) between the upper and lower heating groups in an adjustment range between -50 % and +50 % (not valid for models 30, HPP110eco, IPP110ecoplus, ICPeco, ICP, ICHeco, ICH)						
AtmoCONTROL software ¹¹ for reading out, managing and organising the data logger via Ethernet interface (temporary trial version can be downloaded). USB stick with AtmoCONTROL software available as accessory (on demand)	AtmoCONTROL software ¹⁾ on a USB stick for programming, managing and transferring programmes via Ethernet interface or USB port						
	ControlCOCKPIT with USB port for uploading programmes, reading out protocol logs, activating the User-ID function						
	Displaying of already logged protocol data on the ControlCOCKPIT (max 10,000 values correspond to approx. 1 week)						
Ethernet interface on the rear of the appliance for reading out the protocol log and for online logging	Ethernet interface on the rear of the appliance for reading out the protocol log and for uploading programmes and for online logging						
Double overtemperature protection: Electronic temperature monitoring with freely adjustable monitoring temperature, for models U, I, S with option A6 TWW/TWB (protection class 3.1 or 2), mechanical temperature limiter TB acc. to DIN 12880	Multiple overtemperature protection: Electronic temperature monitoring TWW/TWB (protection class 3.1 or 2 resp. 3.3 for units with active cooling) and mechanical temperature limiter TB (protection class 1) acc. to DIN 12880, AutoSAFETY automatically adjusts to the set value within a freely adjustable tolerance range. Setting individual MIN / MAX values for over/undertemperature and also for all other parameters such as relative humidity, CO ₂						

PID microprocessor control with integrated auto-diagnostic system

Structured stainless steel housing, scratch-resistant, robust and durable; rear of zinc-plated steel

High-temperature connectors on the rear of the appliance for single-phase power connection according to country specific systems and IEC standards

Internal data logger with a storage capacity of at least 10 years

German, English, French, Spanish, Polish, Czech, Hungarian language settings available on the ControlCOCKPIT

Digital backwards counter with target time setting, adjustable from 1 minute to 99 days

The SetpointWAIT function guarantees that the process time does not start until the set temperature is reached at all measuring points – optional for temperature values recorded by the freely positionable Pt100 sensors inside the chamber

> Adjustment of three calibration values for temperature and additional appliance specific parameters directly at the ControlCOCKPIT

As a manufacturer, Menment GmbH + Co. KG clearly labels its devices, which are nectical devices in the sense of the European legislation. The AtmoCOVTROL software is not a medical device.
 All Memment medical devices can be used for their purpose without the software AtmoCOVTROL AtmoCOVTROL is only intended for reading the data logging in conjunction with Memment GmbH + CO. KG medical devices.

SOFTWARE AtmoCONTROL


AtmoCONTROL

The innovative control and logging software

Parameters such as temperature and humidity as well as the process time can be set directly at the ControlCOCKPIT. Ramp programming is done via the control and logging software AtmoCONTROL.

Drag, drop & go!

Numerical and graphic programming of complex processes is a thing of the past. Today, programming is done via AtmoCONTROL by means of the mouse or touchpad on your notebook. Even the most complex ramp programmes are created within minutes. Simply drag & drop the graphical symbols for the desired parameters to the input field and change the values according to your wishes with a mouse click.



Programme functions for appliances with SingleDISPLAY and TwinDISPLAY

- · Reading out, managing and organising the data logger
- · Saving the log memory in various formats
- · Online monitoring of up to 32 connected appliances
- Optical alarms when the alarm limits individually set at the ControlCOCKPIT are exceeded
- Automatic alarm to one or several e-mail addresses

Additional functions for appliances with TwinDISPLAY

- Intuitive programming and archiving of ramps and programme sequences
- Synchronous visualisation of the created programme sequence during programming
- Application-specific repeat functions (loops) can be inserted within a temperature control programme in any place
- Simple creation of repeating weekly programmes
- · Programming, managing and transferring programmes via Ethernet interface or USB port



Device Modifications - Proven and Good

The perfect extension for your Memmert appliance

Our mission at Memmert is to provide you with the best possible solution for your individual application. With the increasing complexity of customer processes, a custom-fit modification of our appliances has many advantages for your application. Through modifications, process and set-up times can be significantly reduced or errors in the application can be completely ruled out by monitoring devices. Even small measures, such as individually adapted accessories, have a noticeable influence on the ergonomics and user-friendliness in the operation of the appliance.

You as a customer have the best ideas - and often already have a specific idea of how our products can be better used in your working environment.

Tell us about your thoughts and let us create an individual solution together with you! Please contact us and call us at +49 9122-925-0 or send us an email to sonderbau@memmert.com.

The Memmert customisation department team is looking forward to hearing from you!

Versatile modifications for our standard appliances



Mechanics

- Customised interior fittings
- Individual entry ports in all sizes and shapes
- . Telescopic slide pull-outs for ergonomic loading



Electronics

- . Extended parameter monitoring e.g. by means of additional measuring sensors
- Optical and acoustic process monitoring e.g. by means of a traffic light system



Software

- Additional interfaces for data evaluation
- Individual temperature, humidity and CO, parameters



Accessories

- Tailor-made subframe and stacking options
- · Modified grids and shelves
- Individual air filters

CUSTOMER SPECIFIC SOLUTIONS

Customised solutions for your requirements

Our expertise as a development partner in plant and project business

The Memmert customisation department has been active in the project business for over 20 years now and has proven itself in countless projects as a strong and reliable partner. The experts in customisation benefit from two aspects: Access to the complete capacities of an ultra-modern and specialised production line, as well as the entire technical know-how of the Memmert company in designing climate and temperature control appliances. Combined with the experience of our project managers, the Memmert customisation department is also able to find a solution for the most complex requirements.

Special sizes

Does your product not fit into a standard unit? We build appliances to measure! Whether you need more volume in the interior or there is not enough space for installation at the installation location, we have the expertise to design your appliance individually. Ask us!

Process and plant integration

Integrate our technology seamlessly into your plant or your work organisation. We will find the right solution together for your process integration:

- · Preparation for integration into your plant
- · Integration of your processes into our appliances
- · Inclusion of customer-specific installations
- . Interface for semi-automatic assembly

Project business

Are you a project developer with ideas for innovative products and looking for a strategic cooperation? Take advantage of our know-how and manufacturing capacities for your project. Our customisation department will be pleased to hear from you!







CLIMATE CHAMBERS

CONSTANT CLIMATE CHAMBERS HPPeto

HUMIDITY CHAMBERS HCP

CLIMATE CHAMBERS (CHero / ICE

ENVIRONMENTAL TEST CHAMBERS CTC / TTO

HEATING AND DRYING OVENS

UNIVERSAL OVENS 1

PASS-THROUGH OVENS UFTS

PARAFFIN OVENS UNpa

VACUUM OVENS VO

INCUBATORS

INCUBATORS I

CO, INCUBATORS ICO

COMPRESSOR-COOLED INCUBATORS ICPeco / ICP

PELTIER-COOLED INCUBATORS IPPeco

MEDICAL DEVICES

JNIVERSAL OVENS Um

NCUBATORS Im

STERILISERS S

CO, INCUBATORS (Comed

BLANKET WARMERS IFEW

WATERBATHS

WATERBATHS WTB

YOUR MEMMERT PARTNER



Memmert GmbH + Co. KG
P.O. Box 1720 | D-91107 Schwabach
Tel. +49 9122 925-0 | Fax +49 9122 14585
E-Mail: sales@memmert.com
facebook.com/memmert.family
instagram.com/memmert.family